## Instructions for use

# ActiGYN<sup>®</sup> VAGINAL GEL

# Moisturizing and soothing vaginal gel.

# PHYSIOLOGICAL pH

# pH 4,5 Gel Formulation

ActiGYN Vaginal moisturizing gel is a medical device with pH 4.5 adjuvant in the prevention and treatment of vaginal and vulvar dryness. With its particular formulation it has a moisturizing, refreshing and soothing effect on the vulvo-vaginal area. Actigyn gel is indicated in cases of vaginal and vulvar dryness and in all situations where an effective soothing effect is needed which gives relief to the discomfort caused by poor intimate lubrication, such as irritation, itching, redness, burning and hypersensitivity or even in case of menopause or if you are taking estrogen-progestins. Hyaluronic acid has a moisturizing and emollient action against vaginal dryness, facilitating the restoration of a natural and pleasant vaginal environment. The lactic acid, a component normally present in the vagina, favours an immediate correction of the pH and facilitates the natural restoration of the delicate vulvo-vaginal ecosystem.

The formulation is enriched by the presence of vitamins with antioxidant action and by plant extracts which, due to their intrinsic properties, exert an emollient and soothing action.

### **HOW TO USE:**

we recommend applying the gel once a day for 6 days or according to medical indications. Before application, wash your hands and affected area thoroughly.

### **FOR EXTERNAL USE:**

Press the tube in order to apply a small amount of gel on the affected area.

### **FOR INTERNAL USE:**

Remove the cannula from the sachet and screw it onto the tube Insert the device thus created into the vaginal cavity

Apply pressure to the tube in order to allow the gel to escape

Unscrew the cannula from the tube

Wash and dry it thoroughly

In severe cases, lubricate the external part before insertion.

### **INGREDIENTS:**

Aqua, Imidazolidinyl Urea, Disodium EDTA, Glycerin, Aloe Barbadensis Leaf Juice, Potassium Sorbate, Ribes Nigrum Fruit Extract, Lavandula Angustifolia Flower Water, Sodium Hyaluronate, Chamomilla Recutita Extract, Lactic Acid, Calendula Officinalis Flower Extract, PEG-40 Hydrogenated Castor Oil, Tocopheryl Acetate (Vitamin E), Polysorbate 20, Linseed acid, Retinyl Palmitate (Vitamin A), Hydroxyethylcellulose.

### **WARNINGS:**


Keep out of the reach of children. Store the product properly sealed and protected from light and heat sources. Do not use the product after the expiry date stated on the package. The expiry date is valid for the correctly stored product. Do not use the product if the package is damaged or open. The device is not sterile. It is not a contraceptive. Do not use in case of organ injury. Do not swallow. Avoid contact with eyes. For pregnant, lactating and young women, the use of the gel must be limited to cases of real need and under the direct supervision of the doctor. Do not use in case of ascertained individual hypersensitivity to the components present. Do not use the product in association with quaternary ammonium compounds. The incorrect observance of the warnings, of the methods of use, as well as the use of the product by several people, could compromise the effectiveness of the product or cause the aggravation of existing problems or cause the occurrence of problems not previously present. Do not use the product for a period longer than 30 continuous days. In presence of an intact hymen, consult your gynecologist before use. It is important to report to the doctor or pharmacist the occurrence of side effects not described in this leaflet.

### **SIDE EFFECTS:**

Under normal and reasonably foreseeable conditions of use, there are currently no known side effects specifically related to the use of the product. Its compatibility has been verified, but in case of irritating or allergic phenomena, stop to use the product and contact your doctor / gynecologist.

### PACKAGING:

30 ml tube.



The medical device must not be disposed in the waste waters or house hold waste. Ask your pharmacist how to dispose the medical devices no longer required. This will help protect the environment.



